





Solutions for EMG /EEG /IOM accessories

# REPUSI CATELOGUE 2025



# ABOUT US

REPUSI Electronics Co., is a manufacturer of Neurology products in China over 20 years. At present ,REPUSI is successfully selling the products in China to large and small hospitals , prestigious medical research institutes ,diagnostic centers and doctors directly or through a nation-wide dealer network.At the same time we expand our service and distribution all over the world .

Looking forward to the cooperation with you ! If any question please feel free to contact us !






# Repusi

REPUSI Disposable concentric EMG needle


Standard sizes , OEM sizes provided.

Disposable Concentric EMG Needle Electrodes	Diameter	Length	Code
	0.35mm (28G)	25mm	Nr. 15.1 Nr. 10.1 Nr. 20.1
	0.45mm (26G)	28mm	
	0.45mm (26G)	38mm	
	0.45mm (26G)	50mm	
	0.50mm (25G)	60mm	

Nr.15.1 serie



Nr.10.1 serie



Nr.20.1 serie






# Repusi

REPUSI Disposable Monopolar EMG needle

Standard sizes , OEM sizes provided.

Disposable Monopolar EMG Needle Electrodes	Diameter	Length	Code
	0.35mm (28G)	25mm	NM3525
	0.45mm (26G)	28mm	NM4528
	0.45mm (26G)	38mm	NM4538
	0.45mm (26G)	50mm	NM4550
	0.50mm (25G)	60mm	NM5060






# Repusi

REPUSI Disposable Injection EMG needle

Standard sizes , OEM sizes provided.

	Diameter	Length	Code	Pack
Disposable Injection EMG Needle Electrodes	0.45mm (26G)	38mm	REPB-38	15 pcs
	0.50mm (25G)	50mm	REPB-50	15 pcs





# Repusi

## REPUSI EMG needle adapt cable

Reusable adapt cable for Disposable concentric EMG needle	Connection plug	Length	Code
	4/5/6-pin	1200mm	GWC-1.2M GWC-1.2F GWC-1.2N






# Repusi

## Disposable Subdermal Needle Electrode

REPUSI offers different range of Disposable Subdermal Needles, with different diameter & length of needle and cable. They are designed to help improve signal quality in electrically hostile environments such as OR and ICU. 12 color leads, 1.5m, 2.0m, 2.5m, 1pc/pack, 24 pcs per box.

	Code	Description
<b>Disposable Sterile Subdermal Needle(single)</b>	REPS1101	Dia 0.4MM *13MM, lead length 1.5M, 12 colors
	REPS1102	Dia 0.4MM *7MM, lead length 1.5M, 12 colors
	REPS1103	Dia 0.35MM *7MM, lead length 1.5M, 12 colors
	REPS1104	Dia 0.4MM *17MM, lead length 1.5M, 12 colors
	REPS1105	Dia 0.45MM *17MM, lead length 1.5M, 12 colors
	REPS1106	Dia 0.40MM *23MM, lead length 1.5M, 12 colors

## REPUSI Subdermal Needle Electrode(single)





# Repusi

## Disposable Subdermal Needle Electrode (moulded head)

REPUSI offers different range of Disposable Subdermal Needles, with different diameter & length of needle and cable. They are designed to help improve signal quality in electrically hostile environments such as OR and ICU. 12 color leads, 1.5m, 2.0m, 2.5m, 1pc/pack, 24 pcs per box.

Disposable Sterile Subdermal Needle(single)	Code	Description
	REPS-03-S2000/0.4*7	Dia 0.4MM *7MM, lead length 2.0M, 12 colors
	REPS-03-S2000/0.4*13	Dia 0.4MM *13MM, lead length 2.0M, 12 colors
	REPS-03-S2000/0.4*19	Dia 0.4MM *19MM, lead length 2.0M, 12 colors
	REPS-03-S2000/0.45*18	Dia 0.45MM *18MM, lead length 2.0M, 12 colors
	REPS-03-S2000/0.45*23	Dia 0.45MM *23MM, lead length 2.0M, 12 colors

## REPUSI Subdermal Needle Electrode(single)






# Repusi

## REPUSI Subdermal Needle Electrode(twist)

### Disposable Twisted Subdermal Needle

Disposable twisted pair subdermal needles are desired to help improve signal quality in electrically hostile environments such as the OR and ICU .12 color leads,1.5m,2.0m,2.5m

Disposable Sterile Twisted Subdermal Needle	Code	Description
	REPS-T1500/0.25/12	Disposable Twisted Subdermal Needle , Dia 0.25MM *12MM length, 1.5M
	REPS-T1500/0.35/7	Disposable Twisted Subdermal Needle , Dia 0.35MM *7MM length, 1.5M
	REPS-T1500/0.40/7	Disposable Twisted Subdermal Needle , Dia 0.40MM *7MM length, 1.5M
	REPS-T1500/0.40/13	Disposable Twisted Subdermal Needle , Dia 0.4MM *13MM length, 1.5M
	REPS-T1500/0.45/17	Disposable Twisted Subdermal Needle , Dia 0.45MM *17MM length, 1.5M





# Repusi

## Disposable Twisted Subdermal Needle(moulded head)

Disposable twisted pair subdermal needles are desired to help improve signal quality in electrically hostile environments such as the OR and ICU .12 color leads,1.5m,2.0m,2.5m

	Code	Description
<b>Disposable Sterile Twisted Subdermal Needle</b>	REPS-03-T2000/0.4/7	Disposable Twisted Subdermal Needle , Dia 0.4MM *7MM length,2.0M
	REPS-03-T2000/0.4/13	Disposable Twisted Subdermal Needle , Dia 0.4MM *13MM length,2.0M
	REPS-03-T2000/0.4/19	Disposable Twisted Subdermal Needle , Dia 0.4MM *19MM length,2.0M
	REPS-03-T2000/0.45/18	Disposable Twisted Subdermal Needle , Dia 0.45MM *18MM length,2.0M

## REPUSI Subdermal Needle Electrode(twist)





# Repusi

## Disposable Single Corkscrew Needle

Disposable single corkscrew needles are desired to pick up EEG or EP signals ..12 color leads,1.5m,2.0m,2.5m

Disposable Sterile single corkscrew Needle	Code	Description
	REPS3301	Disposable single corkscrew Needle , Dia 0.6MM length, 2.0M

## REPUSI Corkscrew Needle Electrode(Single)





# Repusi

## Disposable Twist Corkscrew Needle

Disposable Twist corkscrew needles are desired to pick up EEG or EP signals ..12 color leads,1.5m,2.0m,2.5m

Disposable Sterile twist corkscrew Needle	Code	Description
	REPS4402	Disposable Twisted corkscrew Needle , Dia 0.6MM length, 2.0M

## REPUSI Corkscrew Needle Electrode(Twist)






# Repusi

REPUSI Alligator Clip Wire

## Alligator Clip Lead Wire



- Alligator Clip cable for use with Disposable Surface Electrodes that without leads (see RPSC-S002) .
- Neatly clips onto electrodes, avoiding contact with gel . Availible in red ,green and black
- Ideal for ICU or electrically hostile environment.

	Code	Description
<b>Alligator clip Lead Wire</b>	REP-1.5C	Alligator clips, 1.5M, black/red/green
	GWC/2-1.5C	Alligator clips, cable length 1.2M, extention 20cm (red or black)
	GWC/3-1.5C	Alligator clips, cable length 1.2M, extention 20cm , ( black/red/green)






# Repusi

## Reusable EEG Cup Electrode

Cups are available in gold, Copper and silver/silver chloride ,12 color-coded .

	Code	Description
EEG Cup Electrode	REP-1.5PS	Standard Silver Cup electrode ,Cup diameter :10mm ; Wire Diameter 1.2 mm Length 1.5M
	REP-1.5PC	Standard AgCl Cup electrode ,Cup diameter :10mm ; Wire Diameter 1.2 mm Length 1.5M
	REP-1.5PG	Standard Gold coating Cup electrode ,Cup diameter :10mm ; Wire Diameter 1.2 mm Length 1.5M

REPUSI EEG Cup Electrode






# Repusi

## REPUSI Ground Clamp Electrode

### Ground Clamp Electrode

For use in electrogrgthy as ground electrode , can replace of the Ground strap electrode or Ground plate electrode

	Code	Description
<b>Ground Clamp Electrode</b>	REPG2201	Ground electrode, color red/yellow/green





# Repusi

REPUSI Ground Electrode

## Ground Wrap Electrode

For use in electrography as ground electrode

	Code	Description
<b>Ground Wrap Electrode</b>	REPG3303	Ground electrode (senior), color green, belt 45cm, Snap cable 1.5m





# Repusi

REPUSI Ground Electrode

## EMG Plate Ground Electrode

For use in electrography as ground electrode

	Code	Description
<b>Ground Plate type Electrode</b>	REPG4401	Plate Ground Electrode, cable 1.5m





# Finger Ring Electrode


	Code	Description
Finger Ring Electrode	GWC-1.5R	Finger loops, 5 PIN DIN
	REP-1.5R	1.5mm touch proof socket , 1.5m wire , Paired Ring Electrode





### Disposable surface electrode

The electrodes for recording, stimulation and ground use in EMG/NCS/EP/IONM area .



	Code	Description
<b>Disposable Surface Electrode</b>	RPSC-S002	EKG electrode, 23 x 34mm
	RPSC-S003	Wire attached electrode, four discs, 22mm/25mm .Idea for IOM, SEP, NCS
	NCS-39	Disposable adhesive electrode with lead, 40x40mm
	RPSC-S004	Disposable ground electrode with lead, 60x40mm



# Repusi

REPUSI Handheld simulating Electrode

## Hand held stimulating Electrode

Bipolar hand-held stimulating Electrode ,for NCV studies as stimulating , Electrode center distance:25mm.Round stainless steel electrodes.

Hand held stimulating Electrode	Code	Description
	GWC-1.5S	DIN 5-pin 240° connector,150cm length cable
	REP-1.5S	1.5mm touch proof sockets,150cm length cable
	GWC-1.5S-03	DIN 5-pin 240° connector, 150cm length cable





## Handheld EMG Bipolar stimulator-Width adjustable-25mm / 10mm

For NCV studies as stimulating , 25mm /10MM adjustable from the center to center , 150 CM length cable ,Standard 5 PIN DIN .

	Code	Description
Adjusted Hand-held Stimulator	GWC-1.5S-02	Adult /pediatric
	GWC-1.5H	Adult



## Reusable Recording Bar Electrode

Normally used in nerve conduction and as stimulating points in evoked potential studies.



	Code	Description
<b>Bar Electrode</b>	REP-1.5S-02	Stainless steel $\Phi$ 9mm contacts, 30 mm spacing , 150CM lead cable
	GWC-1.5S-04	20mm Spacing, 150CM cable, Child
	GWC-1.8S-06	40MM Spacing, 150cm cable , Adult





## EEG Ear Electrode ,in pair

Reusable EEG Ear Clips Electrode  
 Standard touch-proof connector  
 Package in pair

	Code	Description
<b>EEG Ear Electrode</b>	REP-1.5E	Cable:1.5M, Material :TPU Connector:DIN 1.5mm Cup:Gold plated copper, 10mm



## Disc Electrodes for NCS

	Code	Description
<b>Disc Electrodes</b>	GWC/2-1.2PC	AlCl1 cups, For Nerve Conduction Study
	GWC/2-1.2PG	Gold-plated cups, For Nerve Conduction Study





### Disposable Ball Tip Stimulation Probe

For intraoperative use ,Strait ,Ball tip,Sterilized with ETO

	Code	Description
<b>Disposable Ball tip Stimulation Probe</b>	REP-1.5N-04	Probe shaftlength:88MM Shaft diameter:1.2mm Tip Diameter :2.0mm Wire length :2m Qty per pack : 1pc

### Disposable Double Ball Tip Stimulation Probes

For intraoperative use ,Strait ,Double ball tips,Sterilized with ETO

	Code	Description
<b>Disposable Double Ball tip Stimulation Probe</b>	REP-1.5N-03	Probe shaft length:80MM Probe Shaft space :5 mm Shaft diameter:1.2mm Tip Diameter :2.0mm Wire length:2m Qty per pack : 1pc





# Repusi

## REPUSI IONM Probes

### Disposable Fork Stimulation Probes

For intraoperative use ,Strait tips,Sterilized with ETO


	Code	Description
<b>Disposable ForkStimulation Probe</b>	REP-1.5N-05	Probe shaft length:90MM
		Shaft diameter:0.8mm
		Wire length:2m
		Qty per pack : 1pc



### Disposable Double Hook Stimulation Probes 180°

For intraoperative use ,Hook ,Sterilized with ETO

	Code	Description
<b>Disposable Double Hook Stimulation Probe 180°</b>	REP-1.5N-06	Probe shaft length:80MM
		Angle :180 degree
		Probe Shaft space :5 mm
		Shaft diameter:1.2mm
		Wire length:2m
		Qty per pack : 1pc






# Repusi

## Disposable Double Hook Stimulation Probes 100°

For intraoperative use ,Hook ,Sterilized with ETO

	Code	Description
<b>Disposable Double Hook Stimulation Probe 100°</b>	REP-1.5N-07	Probe shaft length:80MM Angle :100 degree Probe Shaft space :5 mm Shaft diameter:1.2mm Wire length:2m Qty per pack : 1pc




## REPUSI IONM Probes

## Disposable Double Hook Stimulation Probes ( 90° )

For intraoperative use ,Hook ,Sterilized with ETO

	Code	Description
<b>Disposable Double Hook Stimulation Probe( 90° )</b>	REP-1.5N-08	Probe shaft length:80MM Angle :90 degree Probe Shaft space :5 mm Shaft diameter:1.2mm Wire length:2m Qty per pack : 1pc





### Disposable Flat Tip Stimulating Probes

For intraoperative use ,Hook ,Sterilized with ETO

	Code	Description
<b>Disposable Flat tip Stimulation Probe</b>	REP-1.5N-09	Probe shaft length:90MM Shaft diameter:0.8mm Wire length:2m Qty per pack : 1pc



### Disposable Concentric Stimulating Probes

For intraoperative use ,Hook ,Sterilized with ETO

	Code	Description
<b>Disposable Concentric Stimulating Probe</b>	REP-1.5N-10	Probe shaft length:90MM Shaft diameter:1.3mm Wire length:2m Qty per pack : 1pc





## EMG Laryngeal Electrode

EMG Endotracheal Tube is mainly used to connect with suitable neuromonitor, provide smooth airway for patients and used as medical device to monitor intralaryngeal muscle and nerve activity during operation.



Type	Code	Spec.	I.D.(mm)	O.D.(mm)	Pack
EMG Endotracheal Tube	REPT1201	7.0	7.0	11.0	1pc
EMG Endotracheal Tube	REPT1202	7.5	7.5	11.8	1pc
EMG Endotracheal Tube	REPT1203	6.0	6.0	9.3	1pc
EMG Endotracheal Tube	REPT1204	6.5	6.5	10.0	1pc
EMG Endotracheal Tube	REPT1205	8.0	8.0	12.6	1pc



# Repusi

## Reusable EEG Cap



### EEG accessories

	Code	Type	Head circumference	Age
Silicone EEG Cap	1	N	38-42cm	1-2 new born
	2	I	42-46cm	3-4 Infant
	3	C	46-50cm	5-6 child
	4	S	50-54cm	Small adult, 7-8 child
	5	M	54-58cm	Adult
	6	L	58-62cm	Large adult

	Code	Size
Silicone tube EEG Net cap	REPG-1000	S / M / L




# Repusi

## EEG Bridge Electrodes


	Code	Pack
Bridge electrodes	REP-EEG-101	20 pcs/Bag





苏州乐普斯电子有限公司







Add: Floor 5th, Boji Technology Park ,No.2  
Taishan Road , Suzhou ,China  
TEL:86-0512-68259159 [www.repusi.com](http://www.repusi.com)  
[www.concentricneedleelectrode.com](http://www.concentricneedleelectrode.com)  
[jane@repusi.com](mailto:jane@repusi.com)

